

Voice CPR



Introduction:

This device may Help Save A Life of an employee, co-worker, family member or anyone in the community.

This device is used to assist anyone who is suffering from a Cardiac Event or a Choking situation. It covers Adults, Children and Infants. This device is fully self contained and operates on batteries and needs no internet, cellular or blue tooth to use. Always call 9-1-1 for help First when you begin use of this device to get qualified medical help on the way.

This device is easy to follow by simply pressing a button for the emergency being experienced and follow the human voice instructing you how to help the person.

****This VoiceCPR device by NO means is to replace 9-1-1. It is NOT intended to replace qualified medical help.****

User Instructions:

Read this Manual fully after receiving the device.

When a person is having a Cardiac event or is Choking, Turn the device power switch ON. (It's located on the left side of the case). Then take device to the area where the person is experiencing the emergency and press the appropriate button for the type of medical issue you are dealing with and follow the step by step instructions.

You may Pause and continue play by pressing the Pause/Play button. You can also press the repeat button if you miss or don't understand an instruction. Every press of the Repeat button backs up the instructions by 5 seconds. You can repeat the instructions from the beginning by pressing the first button you pressed.

Continue the instruction procedures until Qualified Medical help arrives.



Battery Replacement:

This VoiceCPR device operates on (4x) AA Industrial Batteries.

The batteries should be replaced annually to assure that it is ready to work at all times.

To replace the batteries, turn the device over and slide the battery door open, remove ALL 4 batteries (always replace all 4 batteries at the same time) and replace with either (4x) AA Industrial Batteries or Class A Alkaline Batteries.

"NEVER USE RECHARGABLE BATTERIES".



Maintenance and Care:

Keep this device in a safe dry location where it is readily available. If an AED is available, place the device on the wall next to the AED. Place it in the optional wall holder or in an easily accessible known location. Check the device monthly by turning ON switch, pressing one of the buttons, then turn OFF.

Gently wipe the device off to remove any foreign debris and keep it clean and ready to use. DO NOT USE ABRASIVE CLEANERS.

Warranty:

This device is warranted for 2 years.

This device will be repaired or replaced due to manufacturing defects or workmanship.

This device is water resistant, but Not Water Proof. This device is not covered if evidence of water damage exists.

Do Not open or tamper with this device or warrantee will be voided.

If failure is experienced, call Hands On The Heart customer service at (727) 433-4542 or Email: HandsOnTheHeart20@gmail.com

Hands On The Heart, Inc.
1423 Wisconsin Avenue
Palm Harbor, FL. 34683

